

SPEC Short Manual

Randolf Butzbach

as of

December 4, 2003

Abstract

This document is intended to provide users with a brief overview of the most commonly used commands in **spec**. **spec** provides users with a powerful language to control their experiment. To program in **spec** is easy and fun. First, the most common macros for positioning, counting, scanning and plotting are briefly described on one page. After a brief introduction in command line editing, a list of commonly used macros follows in a more detailed description. Nevertheless, this list is far from being complete. Each beamline has further standard macros for the control of a particular experiment, which are described in a special section. Then a few notes are given to modify the standard macros, to write own macros and ‘batches’, and how to run them. Even though, it has nothing to do with (but is assumed to be due to) **spec**, a list of frequently observed problems of beamline operation and their solutions is given then. Finally, a list of further help resources is given.

Contents

1	The Most Important Macros	4
2	How to use the Command Line?	5
3	Positioning Macros	6
4	Counting Macros	8
5	Scan Macros	9
6	Plot Macros	11
7	Monochromator Macros	12
7.1	Traditional Style	12
7.2	Pseudo-Motor Style	12
8	Multichannel Analyzer Macros	13
8.1	Tutorial	13
8.2	List of macros	14
9	UNIX/File Macros	20
10	Miscellaneous Macros	21
11	Beamline-Related Macros	23
11.1	Absorption Beamline	23
11.2	Fluo-Topo Beamline	23
12	Writing and Modifying Macros	24
12.1	Macro Commands	24
12.2	How to write my own macros?	26
12.3	How to modify standard macros?	26
12.4	How to write batch files?	27
13	Where to find further macros?	28
14	Trouble Shooting	29
14.1	spec has crashed	29
14.1.1	My motor does not move.	29
14.1.2	When I press the Enter key, I just get blank lines.	29
14.1.3	I cannot abort a movement or scan	30
14.1.4	RST Gamma error: Can't setup connection.	30
14.2	I cannot start spec anymore	31

<i>CONTENTS</i>	3
14.2.1 You have already a valid spec session running	31
14.2.2 You simply closed your spec window	31
14.2.3 You have lost remote control of spec	32
14.3 Further problems when starting spec	32
14.3.1 <code>bind()</code> : address already in use	32
14.4 Workarounds to known bugs	33
15 Further Help-Resources	34

1 The Most Important Macros

All movements and countings can be aborted by typing CTRL-C / STRG-C

shopen [*shutter*]

Opens shutter

shclose [*shutter*]

Closes shutter

umv *motor dest*

Moves motor *motor* to the absolute *dest* (update move).

umvr *motor delta*

Moves motor *motor* relative to current position. (update move relative).

wm *motor [motor2 motor3 ...]*

Prints user and dial position of one or several motors (where motor)

wa

Prints user and dial position of all motors (where all).

set *motor value*

Sets the current user coordinates position of *motor* to *value*

ct [*time*]

Start counting on all counters for *time* seconds. (count)

ascan *motor start finish interv time*

Scans in absolute coordinates.

dscan *motor start finish interv time*

Scans relative to current position.

newfile *filename*

Sets a new data file

cpsetup

Menu for plot parameters

plotselect [*detector ...*]

Selects detector(s) to be plotted.

pplot [*scan_nr*]

Print plot

2 How to use the Command Line?

At ANKA, **spec** is linked with the powerful GNU Readline library, which makes it easy to work with the command line. Here are just a couple of goodies:

- All inputs are saved in a history. You can repeat and edit commands from the history. To see the history, type **history** or **hi** from the **spec** prompt.
- CURSORUP/CURSORDOWN scrolls in the history list.
- PGUP/PGDN + *string* Scrolls in the history, but filters for *string*. *E.g.* when you type `ascan` and then press repeatedly the PGUP-key, you will see all your **ascans** from the past. Clearly, when you make an **ascan**, the change some parameters, issue **ct**, **wm**, *etc* , and then you want to repeat the last **ascan**, simply type `as+PGUP` and you are back to your last **ascan**. You want to have the last but one **ascan**? Press one more time PGUP.
- The TABULATOR-key expands your abbreviated input to the full length. Works for filenames and for the arguments of many (internal) **spec** commands.
- CTRL-R *substring* CTRL-A: After pressing CTRL-R Readline searches online for your *substring* in the whole history. If you want to edit this line before re-launching this command, press CTRL-A first.
- The exclamation mark is substituted by
 - **!!** – the previous command
 - **!*string*** – The most recent command starting with *string*.
 - **!*linenr*** – line number *linenr* in the history.
 - **!*?string*** – the most recent command containing *string*
 - **!**\$**** – the last word in the previous command
 - **!**^**** – the first word in the previous command

There are many more features in readline. Check out the **spec** manual with `h readline` from the **spec** prompt.

Note: You can use all of these tricks on your Linux command line, too!

3 Positioning Macros

mv *motor dest*

Moves motor *motor* (absolute) to the *dest*

mvr *motor delta*

Moves motor *motor* by *delta* mm or degrees

umv *motor dest*

Like **mv**, but the current position is live shown. Useful when moving should be aborted by CTRL-C.

umvr *motor delta*

Like **mvr**, but the current position is live shown. Useful when moving should be aborted by CTRL-C

wm *motor [motor2 motor3 ...]*

Prints user and dial position of one or several motors incl. upper and lower limit

wa

Prints user and dial position of all motors (where all).

wu

Same as **wa**, but prints only user positions of all motors.

ws11 [*all*]

Where slit 1. Macro to display the positions of slit (pseudo) motors in a clear style. Similar macros are available for any other slit (**ws12**, **ws13**, etc). If the keyword `all` is given, also limits and dial values are shown. Actually, any argument will do the trick.

tw *mot [mot2 ...] delta [delta2 ...] [count_time]*

tweaks motors around the current position by *delta* by typing ENTER.

twct *motor stepwidth [count_time]*

tweaks the motor around the current position by *stepwidth* by typing ENTER while showing the current count rate accumulated in *count_time*. When count time is omitted, the default time is assumed.

joy *motor stepsize*

moves the motor joystick-like using the K and L key.

lm [*motor motor2 ...*]

Shows the limits of *motor*. If *motors* is omitted, the limits of all motors are shown.

set_lm *motor low_limit high_limit ...*

The **set_lm** macro is used to establish the low limit and high limit for one motor in units of the user positions.

set *motor value*

Sets the current user coordinates position of *motor* to *value*

set_dial *motor value*

Overwrites the dial (motor controller) position. After a **set_dial**, the absolute position of the axis is usually lost, incl limits! What you want, is probably a **set**, which sets only the user coordinates of your axis, keeping all hard limits in mind.

home *motor [+|-] [home_pos]*

Executes homing on motor *motor*. You can see the current position (according to the controller) by **umw**

Example:

```
home s111; uwm s111
```

Executes homing on slit 1 left blade and shows updated controller register

CEN

Variable, which contains the center of the FWHM of the last scan; can be used in combination with **mv**:

Example:

```
mv th CEN
```

Moves *th* to the center between the upper and the lower position of the half maximum value of the last scan.

COM

Same as **CEN**, but contains the center of mass of the last scan.

usersave

Saves the current motor positions in a file called `recover.mac` and defines a macro **recover** to recover them later.

4 Counting Macros

ct [*time*]

Start counting on all counters for *time* seconds. If *time* is omitted, counting is started for one second. If *time* is negative, counting to $|time|$ monitor counts is enabled. *I.e.* the counters will count until the monitor has seen $|time|$ counts.

uct [*time*]

As **ct**, but the elapsed time and accumulated counts are displayed life.¹

¹Not possible at beamlines where the counter is controlled via Gamma.

5 Scan Macros

There are many different scan commands around for almost all purposes. Here, only the very standard scan commands are listed. All scans can be aborted by typing CTRL-C. When time is negative, the counts are acquired until the detector declared as `monitor` has acquired as many counts as given by `|time|`.

ascan *motor start finish interv time*

scans motor *motor* from *start* deg/mm to *finish* deg/mm (absolute coordinates) in *interv* intervals. Each point is counted for *time* seconds.

dscan *motor start finish interv time*

scans motor *motor* from *start* deg/mm to *finish* deg/mm relative to the current position in *interv* intervals. Each point is counted for *time* seconds.

a2scan *motor1 start1 finish1 motor2 start2 finish2 interv time*

`a2scan` scans two motors, as specified by *motor1* and *motor2*. Each motor moves the same number of intervals with starting and ending positions given by *start1* and *end1*, *start2* and *end2*, respectively. The step size for each motor is $(start-end)/intervals$. The number of data points collected will be $intervals+1$. Count time is given by *time*, which if positive, specifies seconds and if negative, specifies monitor counts.

d2scan *motor1 start1 finish1 motor2 start2 finish2 interv time*

`d2scan` scans two motors, as specified by *motor1* and *motor2*. Each motor moves the same number of intervals. If each motor is at a position *X* before the scan begins, it will be scanned from $X+start$ to $X+end$. The step size for each motor is $(start-end)/intervals$. The number of data points collected will be $intervals+1$. Count time is given by *time*, which if positive, specifies seconds and if negative, specifies monitor counts.

Upon termination, the motors are returned to their starting positions.

a3scan, a4scan a5scan

Analogue **a2scan**, but scans three, four, five axes, respectively.

d3scan, d4scan

Analogue **d2scan**, but scans three and four axes, respectively.

mesh *motor1 start1 end1 intervals1 motor2 start2 end2 intervals2 time*

The mesh scan traces out a grid using *motor1* and *motor2*. The first motor scans from *start1* to *end1* using the specified number of intervals. The second motor similarly scans from *start2* to *end2*. Each point is counted for for *time* seconds (or monitor counts).

The scan of *motor1* is done at each point scanned by *motor2*. That is, the first motor scan is nested within the second motor scan.

A mesh scan creates only one scan entry in the spec data file with a total of $(\text{intervals1}+1)*(\text{intervals2}+1)$ points.

timescan [*counting_time* [*sleep_time*]]

starts the counters in *sleep_time* intervals for *counting_time* seconds. If *counting_time* is omitted, the default time of 1 second is assumed.

rscan *motor start1 end1 interv1 [start2 end2 interv2 ...] time*

allows users to define various measurement density regions among the scanned area; each region is assigned its particular size and intervals number.

lup *motor start finish interv time*

Same as **dscan**, but executes a **mv** *motor* **CEN** after the scan instead of moving to the start position (line up).

6 Plot Macros

newfile *filename*

Sets a new data file

newsample *description*

Defines a new sample description

setplot

Prompts user for plot parameters like on-line update, lin-/log-plot, error bars, etc *for the live plot*.

plotselect [*detector ...*]

Selects detector(s) to be plotted. If *detector* is omitted, the user is prompted for her choice.

Example:

```
plotselect ion1 ion2 ion3
```

Plots first second and third ionization chamber.

cpsetup

Menu to set up the parameters (like titles, data files..) for C-Plot, used by **cplot** and for printing using **pplot**.

cplot [*scan_nr*]

Plot a graph from a datafile on the screen. If *scan_nr* is omitted, the last scan is used.

pplot [*scan_nr*]

Plot a graph from a datafile on the printer. If *scan_nr* is omitted, the last scan is used.

plot3d, contour, pplot3d, pcontour

Analogue to **cplot** and **pplot**, respectively, but makes a 3D or contour plot from a data file.

7 Monochromator Macros

7.1 Traditional Style

Traditionally, the monochromator is driven by a special set of macros:

getE

Returns currently set photon energy.

setE *energy*

Defines current photon energy in keV.

moveE *energy*

Moves monochromator to *energy* in keV.

Escan *start end interv time*

Executes an energy scan from *start* to *end* energy in *interv* intervals.

setmono [*d-spacing*] [*offset*]

sets the *d*-spacing and the vertical beam offset of the DCM. If invoked without arguments, the user is prompted for the values. If the *d*-spacing is given as zero, the user is asked for the crystal type and the Miller-indices. The settings can be saved by using the **sav mono**. Might be disabled to normal users.

7.2 Pseudo-Motor Style

Energy and wavelength can be used like a physical motor, *i.e.* using the usual **mv**, **mvr**, **umv**, **set**, **lm**, **ascan**, ... macros.

Currently, only the pseudo motor **E** is available, which represents a fixed exit setting of the monochromator in photon energy. The unit is keV.

Example:

```
mv E 12.4
```

Moves the monochromator to 12.4 keV (1 Å).

8 Multichannel Analyzer Macros

8.1 Tutorial

1. If not yet running, start **spec**.
2. Start **PyMca &**, a GUI for the MCA, from your Linux prompt.

Attention about the usage of the Graphical User Interface: it is at the moment only capable of monitoring the data acquired, as there is no connection between **spec** and that GUI, appart from the spectrum data itself and its calibration parameters, both stored in shared memory. For instance, the ROIs defined from **spec** or from the GUI are not connected together. They are two separate sets of ROI definitions. The GUI is not capable of sending control commands to the hardware. Be more particularly aware that the acquisition times are not known by the GUI, so if you save the spectrum from the GUI, they wont be included in the file. rather use **spec** for saving your data. The GUI is a display.

3. Run **macsetup** and select the desired parameters.
4. Define the desired ROIs using **mcaroi**. The names of the ROIs are used as reference marks for the definition of pseudo counters.

ROI names	Channels		calibration		ROI	Roi
	first	last	min	max	ACTIVE	numbers
0	0	8191	0.0000	8191.0000	X	0
one	0	1000	0.0000	1000.0000	-	1
two	100	160	100.0000	160.0000	-	2
all	0	8191	0.0000	8191.0000	-	3

5. Start the configuration editor (**edconf**). Type **S** (upper case S) to invoke the scalers page. Insert the defined ROIs as pseudo counters (here: one, two, all):

Scaler (Counter) Configuration

Number	Name	Mnemonic	<>Device	Unit	Chan	<>Use As	Scale	Factor
0	Seconds	sec	ANKA_SC		3	timebase		1e+06
1	Monitor	mon	ANKA_SC	0	0	counter		1
2	Detector	det	ANKA_SC	0	1	counter		1
3	e-Current	ecurr	NONE	0	4	counter		1
4	One	one	NONE	0	5	counter		1
5	Two	two	NONE	0	6	counter		1
6	ALL	all	NONE	1	0	counter		1

6. Acquire data with **ct** and/or any scan command (e.g. **ascan**).
7. Alternatively use **mcaacq [time]** to acquire data. Depending on your settings in **mcasetup**, the data are saved automatically or you have to save your data manually using **mcasave**.

8.2 List of macros

This subsection is taken directly from the original ESRF documentation.

mcsetup [*adc*] [*tmode*] [*gsz*] [*g*] [*cal*] [*Xlbl*] [*Xun*] [*nclr*] [*sleep*] [*bgnd*] [*gui*] [*log*] [*dots*] [*lines*] [*ebars*] [*flag*] [*red*] [*synchr*]

Sets the MCA hardware and software up. You can either give all the parameters on the command line or use the menu that shows up when typing **mcsetup** alone.

Input arguments description. [variable name, default value]

adc: active adc number (1 or 2). [MCA_ADCNO, 1]

tmode: 1 if live time, 2 if real time. [MCA_TMODE, 2]

gsz: MCA memory group size in channels. [MCA_MEMGRPSIZE, 8192]

g: MCA selected memory group. [MCA_MEMGRP[1], 0]

cal: 1 if spectrum is calibrated, 0 otherwise. [MCA_ENERGY, 0]

Xlbl: X plot label for calibrated spectrum. [MCA_CP_XLBL, ``Energy``]

Xun: X unit for calibrated spectrum. [MCA_CP_XUNT, ``KeV``]

nclr: 0 for memory clearing prior to acquisition, 1 otherwise. [MCA_SAVEBUFFER, 0]

sleep: number of sec. sleeping between run plot updates. [MCA_DISP_SEC, 0]

bgnd: 1 if you want background subtraction in ROIs counts. [MCA_BACKSUB, 0]

gui: 1 for gui plot. [MCA_GUI, 0]

log: 1 for log plot. [MCA_LOG, 0]

dots: 1 for large dots in plot. [MCA_DOTS, 0]

lines: 1 for lines plotting. [MCA_LINES, 0]

ebars: 1 for error bars plotting. [MCA_EBARS, 0]

flag: data saving flag [MCA_FLAG, 0]

0x01 save spectrum during scans to a file.

0x02 save spectrum after each other acquisition to a file.

0x04 save to specific MCA file, otherwise use scans file [DATAFILE].

red: data reduction coefficient, to be applied at saving. [MCA_REDUCTION, 1]

synchr: MCA acquisition is hard-synchronized with **spec** configured timer. [MCA_SYNCHRO, 0]

Note: It is possible to have the last acquisition elapsed times values reported in counters by simply configuring the counters you need out of `mcaLt`, `mcaRt` and `mcaDt` mnemonics, respectively standing for Live, Real and Dead times.

mcaoff

Disables the MCA. The MCA is no longer addressed during standard `spec` operation (`ct`, `scans ..`), but still from `mcaacq`.

mcaon

Enables back the MCA.

mcainit

Initialises the MCA macros and hardware. Can be called at startup to re-setup things without giving any parameter, nor calling any menu. The current or default parameters are then set. Default are mentioned together with `mcasetup` description.

mcaacq [*preset_seconds*] [*roi | first*] [*last*]

Acquires data with the MCA only. The rest of counters configured in `spec` are not addressed by this command. If you need it, use standard `ct` command. (don't forgot, `mcasetup` set AUTO-RUN mode ON, or `mcaon`). `mcaacq` consists in the following sequence:

- stopping any already running acquisition.
- clearing (if setup) the active MCA memory group.
- presetting the time.
- starting the acquisition.
- polling the device till time is over.

meanwhile:

- reading the MCA buffer.
- plotting the data spectrum.

Saving: (if set-up) the acquired spectrum to disk file.

The MCA memory is read within the range specified, or if none, the active ROI. If the requested preset time is 0 or un-specified, the acquisition is started "for ever"; press CTRL-C or "s" or "q" to stop it. Pressing "c" would clear the memory on the fly. If the GUI is not active, you can also change the plot update intervals (press "u"), change the plot x or y ranges (press "x" or "y"), integrate counts on region (press "i"), toggles plot attributes such as lines (press "l"), large dots ("d") and logy ("g"). Meanwhile, the acquisition goes on.

mcaclear

Clears the active MCA memory group.

mcastart [*preset_seconds*]

Presets and starts an acquisition.

mcastop

Stops acquisition if running.

mcahalt

Send “stop acquisition” signal whatever is the status of the device. It might happen that the device server gives a meaningless error message, if the status is already “IDLE”.

mcaread [*roi* | *first*] [*last*]

Reads the MCA memory within the range specified.

mcatimes

Reports the MCA elapsed live, real and dead times.

mcastat

Prints out a detailed status of the **spec** MCA registers.

mcaroi [*nr* | *name*] [*first last*] [*name*]

ROIs definition.

detail of usage:

mcaroi nr first last counter

modifies ROI number *nr*.

mcaroi nr

sets ROI *nr* as the active one.

mcaroi nr counter

associates name “counter”, which is eventually a counter mnemonic in config, with ROI number *nr*.

mcaroi

calls the ROI menu.

#@CTIME 1 17 17

Time preset, MCA elapsed live time, MCA elapsed real time.

#@CALIB 0 1 0

Calibration parameters as $a b c$, in $E = a + b \times \text{channel} + c \times \text{channel}^2$.

@A 0 0 0

MCA data. Each value is the content of one channel, or an integrated value over several channels if a reduction was applied.

Data reduction is useful in some cases to minimize file sizes, which might grow very fast and eventually fill up the disk. It consists of:

- averaging counts every *factor* points.
- multiplying by *factor* each integer average to get an integrated value.

When `-silent` the macro does not ask for a user comment.

mcaload *file scan_no [scan_point_no] [0|1]*

Loads spectrum from a file. The loaded spectrum is stored into WS memory as if it would have just been acquired. If file is a standard scan file, spectra are referenced to by their *scan_no*, and eventually their *scan_point_no*. Both numbers must be 0 when loading data from private MCA files. 4th argument to 1 means that the calibration parameters are to be loaded from the file as well. In that case any other current calibration is overwritten. To avoid it, **mcasave** it with the current spectrum, before **mcaload**ing. Reduced data channel numbering starts from the first channel value (prior to reduction) added integer part of half the *factor*, every next point incremented by *factor*.

mcapar

Shows calibration fit parameters or inputs new parameters.

mcaE

Toggles between channel (uncalib) and calibrated mode.

mcacal

Computer aided energy calibration of the MCA.

mcacalm

Manual calibration of the MCA.

mcapeak

Peak search on the current spectrum.

mcashow

Shows peaks found.

mcacplot *[first|roi] [last]*

A screen plot of the current MCA data in cplot format. Part of the plot parameters are taken from **mcasetup**, rest from **cpsetup** settings.

mcapplot *[first|roi] [last]*

A printer plot of the current MCA data in cplot format. That follows the same rules as **mcacplot** macro.

mcaguion

Turns MCA Graphical Interface on.

mcaguioff

Turns MCA Graphical Interface off.

mcasplot *[xmin] [xmax] [ymin] [ymax]*

Plots the current MCA data spectrum. Does nothing when using GUI.
Detail of usage:

mcasplot first last [ymin ymax]

plots acquired spectrum between the limits specified. Eventually, *ymin* and *ymax* are *Y* limits values so that you may *Y*-rescale your plot.

mcasplot roino [ymin ymax]

plots acquired spectrum within ROI specified.

mcasplot

plots the whole acquired spectrum

9 UNIX/File Macros

newfile *filename*

Sets a new data file name

fon *filename*

switches logging of all screen input/output in file *filename* on. **foff** switches logging off again.

foff *filename*

switches logging of all screen input/output in file *filename* on.

cd *directory*

Change directory. Without argument, cd change to the home directory

pwd

Print working directory

ls [*reg_exp*]

list files in working directory

l [*reg_exp*]

Abbreviation for **ls -l**; list files long (detailed) form in working directory

unix("command"), u *command*

sends any command to the UNIX operating system

10 Miscellaneous Macros

shopen [*shutter*]

Tries to open shutter *shutter*. If *shutter* is omitted, the last used shutter is opened. If opening is forbidden by the safety system, a message is shown. *shutter* is one of the following mnemonics (aliases are given in parentheses):

Absorber	:	abs (0)
Front-End Shutter	:	fs1 (1, fs)
Safety-Shutter	:	fs2 (2, zs)

shclose [*shutter*]

Closes shutter; if *shutter* is omitted, the last used shutter is closed. Mnemonics are same as in **shopen**.

shstate

shows the state of all shutters.

init_gamma

Initializes all RST Gamma ANKA Motors on the OS/9 System. Necessary only after a system reset of the VME crate.

restore_gamma

Restores the RST Gamma parameters from Gamma's database. Useful when working with other measurement software which relies on Gamma's parameters: Issue a **restore_gamma** before leaving **spec** and starting the third-party application.

reset_gamma

kills and restarts all notorious problematic processes of Gamma (currently: network and CAN-bus). Requires a responding TCP/IP network on the OS/9 side, which might be disabled from the Gamma modules, too. The communication output is written to `/tmp/reset_gamma.log`

cryforhelp

Sends an e-mail to computing services for assistance.

print, p *expression*

prints *expression* to the screen. Expression can be also a mathematical expression, so that you can use **spec** as a calculator:

Example:

```
74.SPEC> p dhkl=5.431/sqrt(pow(3,2)+pow(1,2)+pow(1,2))
1.63751
```

```
75.SPEC> p deg(asin(1./(2*dhkl))), "deg"
17.7787 deg
```

Calculates the Bragg angle in degrees for the Si(311) reflection and 1 Å wavelength

comment, com *comment*

Writes comment to data file.

sleep *time*

Sleeps for *time* seconds.

config

Invokes the config-editor.

history, hi

shows the last submitted commands.

Commands can be repeated by

!linenr

or

!abbreviation

Examples:

!156 Repeats command number 156

!dsc Repeats last command which starts with **dsc** (e.g. **dscan ...**)

For a full description of what is possible see `h readline`

sync

Synchronizes **spec**'s internal motor data base of with the controller values. If there is a discrepancy, **spec** prompts the user whether the controller (dial) value should be overwritten. If answered with no, **spec** aligns its database to the dial value.

emac *macroname*

Edits the existing macro *macroname* with the default editor (defined in the variable EDITOR) and loads the changed macro in **spec**..

11 Beamline-Related Macros

11.1 Absorption Beamline

kscan *energy_motor edge_energy time*

Executes a scan around the absorption edge at *edge_energy*. The detailed scan can be defined by using the **kscan_setup** macro. Default values can be used using the **kscan_default** macro.

Currently, only **energy_motor** is available: **E** executes a fixed exit scan.

11.2 Fluo-Topo Beamline

fsopen [*time*]

opens the fast-shutter at the topography station for *time* seconds. If *time* is omitted or equal zero, the shutter is open for infinite.

fsclose

closes the fast-shutter at the topography station after being open for infinite.

fsstate

shows the state of the fast-shutter and the safety-shutter

12 Writing and Modifying Macros

Almost all commands you issue to **spec** are macros. The **spec** programming language is very similar to the C-language; so most of the people will feel immediately at home. This section describes the differences to C and things you need to know for programming **spec**. For help, it is a good idea to look in the definition of other macros using **prdef**. For a detailed description of the **spec** language refer to the **spec**-Homepage (www.certif.com)

12.1 Macro Commands

prdef *macro*

Prints the definition of *macro*

def *macroname* ' *statements* '

Defines *statements* as macro.

cdef ("macroname" , *statements* [, *key* [, *flags*]])

Chain definition: Appends *statements* to an (existing) macro *macroname*. With the optional *key* argument, the pieces can be selectively replaced or deleted, *i.e.* by using the *key*, parts can be later accessed. The *flags* argument controls whether the pieces are added to the front or to the back of the macro or whether the pieces should be selectively included in the definition based on whether *key* is a currently configured motor or counter mnemonic. The bit meanings for *flags* are as follows:

0x01 : only include if *key* is a motor mnemonic and the motor is not disabled.

0x02 : only include if *key* is a counter mnemonic and the counter is not disabled.

0x10 : place in the front part of the macro.

0x20 : place in the back part of the macro.

If *flag* is the string **"delete"**, the piece associated with *key* is deleted from the named macro, or if the name is the null string, from all the chained macros. If *key* is the null string, the **FLAGS** have no effect.

For an introduction in using the **cdef** function see the ESRF tutorial (<http://www.esrf.fr/computing/bliss/tutor/spec.html>).

undef *macro*

Undefines *macro*

lsdef [*reg_exp*]

lists all macros currently known to **spec**.

qdo *file*

Executes macro/batch/script from file *file*. Also to be used to load a macro, if *file* contains a macro definition, *i.e.* the definition of a macro is executed.

savmac *macro_name file_name*

Write macro *macro_name* to file *file_name*

jtdo("*macro*")

Loads the macro *macro* from the installed standard macro collection. **jtdo** searches first in the system user directory, then in the beamline directory, then in the general ANKA directory, and finally in the ESRF macro directory.

Note, that we have most of the ESRF macros installed here, too. Many of them can be used out of the box, others contain special ESRF commands (**esrf_io**) or hardware and are of no direct use at ANKA. However, they can be useful as a template or to get an idea for your own macros.

emac *macroname*

Edits the existing macro *macroname* with the default editor (defined in the variable EDITOR) and loads the changed macro in **spec**..

moredef *macroname*

Like **prdef** but uses a pager like more to display the macro definition. You can define PAGER to change the page to less or something else.

#

Starts a comment until the end of the line

\$#

Replaced by the number of arguments given.

\$1

Replaced by the first argument given, when macro is invoked.

\$*

Replaced by all arguments given, when macro is invoked.

12.2 How to write my own macros?

1. Open/create a file with your favourite editor. From the **spec**-prompt, you can use *e.g.* `vi mymacro.mac`. If you are not familiar with `vi`, try something like `nedit`, `emacs`, `xemacs`, `joe`, `pico`, ...
2. Write the macro starting with the keyword **def**, then give the macro name and enclose your code in `' '`, *e.g.*

```
def hello '
    print "Hello World!"
'
```

3. Save the file
4. Load the file in **spec** with **qdo**, *e.g.*

```
1.SPEC> qdo mymacro.mac
```

5. Run the macro from the **spec**-prompt:

```
2.SPEC> hello
Hello World!
```

12.3 How to modify standard macros?

1. Search for the desired position to be changed, using the **prdef** command. Example (to change the update value format when using **umv**):

```
3.SPEC> prdef umv

# SPECD/standard.mac
def umv '_mv $*; _update("$1")'
```

2. Save the standard macro to your own file, *e.g.* using the **savmac** macro: Example:

```
4.SPEC> savmac _update myupdate.mac
```

3. Proceed with Section [12.2](#).

Hints:

- The macro **emac** (edit macro) does all of this on the fly:

Example:

```
emac _update
```

Creates a temporary file, opens an editor (defined in the variable `EDITOR`), and loads the macro afterwards.

- To restore the original **spec** macros, type **newmac**.

12.4 How to write batch files?

The description of how to write macros in Section 12.2 describes actually already how to write a batch file. The command `qdo file` executes everything what is written in that given file: After writing our macro in Section 12.2 we have executed the macro definition. However, we can put anything else instead of or additionally to the macro definition. The file could read for instance like this:

```
def hello '  
    print "Hello World!"  
,  
hello  
ascan th 0 10 20 1  
mv th CEN  
dscan th -1 1 100 3  
pplot
```

First we define the `hello` macro (`def ...`), which we execute immediately after definition (`hello`). After the greeting, we do a coarse absolute scan (`ascan ...`) of `th` from 0 to 10 deg. Next we move `th` to the center of the reflection (`mv th CEN`) and then we continue with a fine scan around the center position at which `th` is now positioned (`dscan ...`). And finally, we send the plot to the printer (`pplot`).

13 Where to find further macros?

`/usr/local/lib/spec.d/local/user` – where you can put your own macros to provide them to other users

`/usr/local/lib/spec.d/local/BL` – Beamline specific macros

`/usr/local/lib/spec.d/local/anka` – ANKA standard macros

`/usr/local/lib/spec.d/local/esrf` – ESRF standard macros

<http://www.esrf.fr/computing/bliss/spec/local>

14 Trouble Shooting

Please note: *spec* has proved to be very stable. Problems that occur are likely not to be due to *spec*, but due to the underlying Gamma control system. Moreover, *spec* is designed to be failsafe. I.e. the idea of simply closing your *spec* window will not help but rather bring you in deeper troubles!

14.1 spec has crashed

Are you sure, that it is really *spec* which crashed? Probably, the underlying control system Gamma has crashed und you just see *spec* waiting for the underlying control system, which does not respond anymore. Check out the list below for detailed help.

14.1.1 My motor does not move.

There are several reasons for this phenomenon:

- **You control the motors via Gamma and have forgotten to initialize the motor driver:**

Enter `init_gamma` from the *spec* prompt. Depending on your beamline, you might need to wait about one minute for the sysload on the OS/9 computer to decrease. If you continue your work too soon, Gamma will crash and you need to reset the whole system.

- **Your Gamma driver for your motor controller has crashed.**

Reset your Gamma Motor Driver as described in Section 14.1.2. Ask your beamline scientist for assistance.

- **Your McLennan PM595 Motor Controller has crashed.**

Reset your PM595 Motor controller. Ask your beamline scientist for assistance.

- **Due to a hardware exception, both limit switches appear to be active and spec does not inform you about this state.**

Check out, whether all cables are mounted and the opto-coupler board is switched on. Check the communication between *spec* and your motor controller (switch on debugging with `debug 192`, switch back to normal operation with `debug 0`).

14.1.2 When I press the Enter key, I just get blank lines.

Probably, your CAN-Bus interface in Gamma has crashed. Try to run `reset_gamma` from your Linux prompt. `reset_gamma` will kill all problematic processes on your

underlying control system and restart them. You also have to run `init_gamma` from your `spec` session afterwards, in order to be able to use your motors again.

If this does not help, ask your beamline scientist for assistance.

14.1.3 I cannot abort a movement or scan

You have typed (several times) CTRL-C, but nothing happened or `spec` replies with

```
Waiting for motors to stop.
Still waiting.
```

or something similar. Probably, your motor driver in the control system or your motor controller (likely at McLennan PM595) itself hangs. In the former case, you need to reset the driver in Gamma as described in Section 14.1.2. In the latter case, you need to reset the motor controller. In both cases, ask your beamline scientist for assistance!

14.1.4 RST Gamma error: Can't setup connection.

You get an error message like:

```
RST Gamma error: Can't setup connection.
RST Gamma unavailable.
```

You have lost your network connection to the Gamma control system. (Perhaps you have closed and restarted your `spec` session without waiting at least 15 seconds.) Try the following:

- **Wait for at least 2–3 minutes** (check with a watch!) – Then try again. Sometimes the network interface becomes available again.
- **Restart the Gamma network interface automatically** (works only, if you still can make a `telnet` connection to the OS/9 computer:
 1. Run `reset_gamma` from the Linux prompt. If this does not work, skip the rest of this item and proceed with “Restart the network interface manually”.
 2. Run `anka_par("reconnect")` from the `spec` prompt.
 3. Run `init_gamma` from the `spec` prompt.
- **Restart the Gamma network interface manually** (You might need assistance from your beamline scientist):
 1. Open a terminal window on the OS/9 computer (click with the mouse in the backdrop) and enter the password (ask your beamline scientist).
 2. Find out the node- and process number (PID) of your network interface (note that the UNIX “|” is a “!” in OS/9):

```
GAMMA: procs ! grep CN
103  0  0.0    128 112.00k  0 s 16:08:03.80 1294:21 CNetMan <>>>nil
106 103  0.0    129  80.00k  0 s 79:41:28.46 1294:21 CNETTCP_309 <>>>nil
```

Here, 103 and 106 are the PIDs and 309 is the node number.

3. Type the following, where you replace the example PIDs and Node numbers by your respective PIDs and Node number:

```
GAMMA: setenv LNODE 309
GAMMA: kill 103 106
GAMMA: cnetman <>>>/nil &
```

4. Run `anka_par("reconnect")` from the `spec` prompt.

14.2 I cannot start spec anymore

You want to start `spec` and you get something like

```
Can't lock state file "/usr/local/lib/spec.d/spec/userfiles/sul_ttyp#L".
Are you already running on this terminal or another virtual tty?
```

There are several causes that you can not start `spec` anymore:

14.2.1 You have already a valid spec session running

Close your running `spec` session first or continue with your running `spec` session: There are several beamline control components which can be used by one process (*i.e.* `spec` session), only. Thus you can not run more than one `spec` session at once which uses the same control components.

14.2.2 You simply closed your spec window

`spec` is programmed to be failsafe. Closing your `spec` window will *only* close the window while `spec` considers this as a failure on the user interface side (like the breakdown of a network connection, too) and continues the operation. In other words: when you lost your network connection, your long-run scan will not be aborted.

Unfortunately, you don't have a regular user interface to `spec` anymore. So you have to tell `spec` to terminate by operating system signal:

1. Lookup the process identifier (PID) of your `spec` session:

```
linux:~> ps x | grep spec
18822 pts/6    S          0:00 spec
```

where you need to replace the string `spec` after `grep` by the name of your `spec` session (*e.g.* `fourc`, `optics`, *etc*).

2. Tell **spec** to terminate ordinary:

The following command will tell **spec** to terminate normally (*i.e.* saving all files *etc*):

```
linux:~> kill -HUP 18822
```

and check, whether **spec** really has terminated, by repeating the `ps` command above. If **spec** has terminated, stop here. If not, try the next point:

3. Tell **spec** to end:

```
linux:~> kill 18822
```

Check again whether **spec** is still running or not (*c.f.* first item). If **spec** is still running, try the last alternative:

4. Tell Linux to kill **spec**:

You should consider this option as the very last way. It is equivalent to switching off your computer. **spec** will not learn about your termination request. Thus, your variables and history *etc* will not be saved!

```
linux:~> kill -9 18822
```

14.2.3 You have lost remote control of spec

You are running **spec** over a network connection (*i.e.* `ssh` or `telnet`) and you have lost/closed this connection. Now, you can not access **spec** anymore:

Log on to the computer on which **spec** is running and proceed with Section [14.2.2](#).

14.3 Further problems when starting spec

14.3.1 `bind()`: address already in use

You run more than one program (**spec** session? Status display?) on the very same computer which talks to the very same OS/9 computer (Gamma control system).

Sorry, that's not possible. Close one of those programs. Ask your beamline scientist for assistance.

14.4 Workarounds to known bugs

Unfortunately, ANKA's control system and hardware suffer from several well-known bugs. ANKA is striving to get these bugs removed. Until then, there are a couple of work-arounds of these bugs.

Connection to beamline: (Applies, when using RST's Gamma for beamline control)

When you close your **spec** session, take care that you *wait at least 15 seconds* before you restart **spec**. Otherwise, the entire network interface of the VME computer will crash and your beamline scientist has to restart the network on location; there is no remote access possible anymore.

Small step discrepancy when moving by small distances: (Applies, when using the McLennan PM595 motor controllers) When you want your motor to move for less than five steps, the McLennan motor controller will return that he is already at this location and will not move. After the movement, **spec** will compare its own target position with the target position reached. Since this discrepancy becomes obvious here, **spec** will ask who is right and who is wrong. Usually, it seems to be a good idea to say no to the question

```
Should the ANKA Gamma Motor controller be changed [spec suggests NO]?
```

However, it was also observed that the motor controller said that there is nothing to move, but the motor moved nevertheless. In this case, you should say yes, in order to update the controller to the physical position of the motor.

15 Further Help-Resources

- **help** *topic*

Online help from the **spec**-command line: Example: `h readline` – shows the help page for the command line editor. When invoked without argument, a list of topics is shown.

- **cryforhelp**

sends e-mail to computing services.

- <http://www.certif.com>

The official **spec** home-page

- <http://www.esrf.fr/computing/bliss>

The ESRF-BLISS Homepage

- <http://www.esrf.fr/computing/bliss/spec/local>

ESRF-**spec** macros; help and downloading

- <http://www.esrf.fr/computing/bliss/tutor/commands/commands.html>

spec commands reference list

- <http://www.esrf.fr/computing/bliss/tutor/spec.html>

Tutorial in **spec**